
DENTSPLY

SmartLite[®]
max
L.E.D. Curing Light

DIRECTIONS FOR USE - ENGLISH
INSTRUCCIONES DE USO - ESPAÑOL
MODE D'EMPLOI - FRANÇAIS

TABLE OF CONTENTS

English..... 1

Español 10

Français..... 20

Caution: For dental use only.
U.S.A.: Rx only.

1. PRODUCT DESCRIPTION

SmartLite® Max L.E.D. Curing Light is a pen-style, high-power output, cordless light polymerization device for use in dentistry. SmartLite® Max Light offers a dual advantage of super fast and broad spectral curing of resins. SmartLite® Max Light incorporates two types of LEDs to achieve broad spectral output to polymerize virtually all the resins in the market and a high power mode to produce super fast curing effect. It is a small probe-style light with many features normally found on larger lights. The curing performance is proven to be similar to the conventional halogen curing light, making it a versatile LED curing light suitable for virtually all chairside curing light tasks.

1.1 Product components

Handpiece Assembly
Curing Probe housing the LEDs
Charging Base Assembly
Wall Plug-in Power Supply
Eye Protection Assembly
Barrier Sleeves

1.2 Product overview

The SmartLite® Max Light utilizes state-of-the-art semiconductor diodes to generate the light energy used by this device. This technology allows for a more reliable product that uses less energy and is housed in a smaller package than traditional halogen lamp curing lights. This allows the creation of a device that is small, portable, efficient, and should provide many years of consistent use. SmartLite® Max Light features a built in rechargeable battery, allowing cord-free operation. The power supply can be attached if needed, ensuring consistent uninterrupted use the entire day. The SmartLite® Max Light also features a fully rotating head, with a complete 360° range of movement, allowing you to see the controls while applying light at any angle desired.

SmartLite® Max L.E.D. Curing Light features:

- Probe housing the 2 LED's
- High power output (1200-1400mW/cm²) ensures super fast 5-second cure for some commonly used (A3 and below) composite resin shades (see Curing Card for specific times)
- Wide spectral output (400-460nm) cures virtually all available VLC materials
- Cordless with corded back up for uninterrupted use
- 360° rotating probe
- Sleek, convenient and lightweight (220g/7.8oz)
- Quiet, no fan
- Rechargeable lithium ion battery
- Simple controls, push buttons labelled with icons
- Convenient, informative, attractive display with icons for different modes, curing time selection and the battery status
- Four cure modes, easy to select, with mode icon on display (standard, boost, ramp and pulsed)
- 90° angle light guide
- Automatically goes to sleep mode in 10 minutes
- Adjustable timer (5, 10, 15, 20 seconds) with countdown display
- Audible signals confirm button presses to select the mode and time, and every 5 seconds during the curing process
- Built in radiometer on base with color-coded LED scale indicator
- The unit turns off if the temperature at the probe tip reaches 41°C to minimize exposure of personnel and patients to excessive temperatures
- Sustained high energy light output regardless of the battery charge during the operation

1.3 Indications

Curing of all dental products with camphorquinone (CQ) and L-TPO photoinitiators.

1.4 Contraindications

Use in patients prone to photobiological reactions: Do not use the SmartLite® Max L.E.D. Curing Light on patients with a history of photobiological reactions (including patients with Urticaria solaris or erythropoietic protoporphyrie) or those currently undergoing treatment with photosensitising pharmaceuticals.

1.5 Technical specifications

PARAMETER	VALUE
Dimension in mm (length x diameter)	240 x 30
Range of Light	377-490nm
Output at the Tip (mW/cm²)	1200 (standard)/1400 (boost)
Weight in Grams	220
Goes to Sleep Mode in	10 minutes
Battery Life	At Standard mode, over 200 cycles of 10 seconds
Time for Complete Charging	< 3 hours
Battery	Rechargeable Lithium Ion
Operating Temperature	Unit turns off automatically if the temperature at the probe tip reaches 41°C
Power Input, Curing Light	12VDC @ 2.0A
Power Input, Power Supply	90-240VAC @ 50-60 Hz
Pulse Rate	10 Hz
Pulse Duty Cycle	50%
Operating Temperature	15-30°C (59-86°F)
	0-95% RH
	-1000 - +8000 meters elevation

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these directions.

SAFETY SYMBOLS

A number of symbols are used both on the device and throughout this manual. It is important that you understand the meaning of each one, so that when encountered you will be familiar with the type of information that is being expressed. Obey all safety messages that follow this symbol to avoid possible injury.

Safety alert symbol
This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

Protect your eyes from the light generated by this product.

Beware of electrical shock or injury.

Type B Applied Part: Protection from electric shock: Type B

2.1 Warnings

1. The SmartLite® Max L.E.D. Curing Light emits high intensity visible light radiation. Prolonged continuous exposure to eyes, oral mucosa or exposed skin can result in injury or irritation. When curing, position the light probe tip directly over the polymerizable material and avoid extraneous exposure to soft tissue for prolonged periods. Wear suitable protective clothing and eye protection designed to filter blue-violet and ultraviolet radiation. Securely attach protective light shield before use.

Eye exposure: Avoid prolonged exposure to unprotected eyes. Prolonged exposure may cause blurring or starring of vision, or result in visual after-image. In most cases these phenomena are temporary and do not require treatment. If vision abnormalities persist, consult a physician.

Skin contact: Avoid prolonged exposure to skin. If skin rash, sensitization or other reaction occurs, discontinue use and seek medical attention.

Oral mucosa contact: Avoid prolonged exposure to oral soft tissues. Prolonged exposure may cause soft tissue irritation or burns. Most minor reactions require only thorough cleansing and palliative treatment. If irritation persists, seek medical attention.

2. The SmartLite® Max L.E.D. Curing Light is an electric device. To avoid electric shock:
 - Do not attempt to change the Curing LED, open or alter the unit in any way (other than as directed in these Directions for Use).
 - Do not introduce any foreign objects into the housing of the unit.
 - Connect the power plug into a suitably grounded and approved outlet. If using an extension cable, make sure that the ground line is not interrupted.
 - Always unplug the base/charger before disinfecting.

3. The SmartLite® Max L.E.D. Curing Light is not suitable for use in the presence of a flammable anesthetic mixture with air or nitrous oxide. Use in well ventilated area.

2.2 Precautions

1. This product is intended to be used only as specifically outlined in these Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and is the sole responsibility of the practitioner.
2. DO NOT look directly at the light emitted from the curing probe. As with any bright light source, staring directly at the source may cause temporary effects such as after-image. (See Warnings)
3. It is recommended that all curing light users and patients employ suitable eye protection during operation of the SmartLite® Max Light. (See Warnings)
4. Persons who have had surgery for eye cataracts are especially sensitive to light and generally should not use or be exposed to the SmartLite® Max Light unless they employ suitable protective measures. These measures should include the use of spectacles that filter out blue-violet and ultraviolet light. (See Warnings)
5. Safety of the SmartLite® Max Light has not been established in people who have been fitted with an implanted heart pacemaker and have been cautioned against the use of small electrical appliances (such as electric shavers, hair dryers, etc.). In these cases we recommend that the SmartLite® Max Light not be used.
6. DENTSPLY provides light probes suitable for use with the SmartLite® Max Light. Other light probes must not be used.
7. Curing probes must be maintained in a clean, undamaged condition. Probes that are scratched or damaged in any way should not be used.
8. Use of a protective barrier sleeve over the probe lens is strongly recommended. Do not allow contact of restorative materials with probe lens. Attempts to remove cured material from lens may scratch lens. A contaminated or chipped/scratched lens will reduce the effectiveness of the unit.
9. Devices marked "Single Use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
10. Check unit for adequate light output before each procedure. Failure to verify output may allow inadequate curing (see Radiometer Operation section). Only the built-in radiometer should be used. Other radiometers may give erroneous output data.
11. DO NOT IMMERSE UNIT IN WATER OR DISINFECTANT. Do not spray liquids directly onto the light. Spray a towel, then wipe the light. Prevent liquids from entering openings on unit, especially the socket for the probe. Always wipe the unit with the probe installed.
12. Should the SmartLite® Max Light become damaged in any way, do not attempt to operate the unit unless thoroughly inspected and serviced by qualified DENTSPLY Service personnel.

13. The SmartLite® Max Light does generate some heat when in use. If multiple curing cycles in rapid succession occur, especially in Boost mode, the case may get warm. The unit is equipped with an automatic thermal safety shut-off and will shut down if the case temperature reaches 41°C. If this occurs, to avoid skin or mucosa burns, place the unit in the charging cradle or on the counter and wait at least 6 minutes before trying to use the unit again.

2.3 Adverse reactions

1. Prolonged unfiltered exposure to the light source may cause damage to the eye. (See Warnings)
2. Prolonged contact with soft tissue may cause injury or irritation to the tissue. (See Warnings)
3. Medical conditions such as solar urticaria, erythropoietic protoporphyria or cataract surgery may be aggravated by exposure to emitted light. (See Contraindications, Precautions)
4. Electrical appliances can interfere with some cardiac pacemakers. (See Precautions).

2.4 Storage conditions

Inadequate storage conditions may shorten the useful life and may lead to malfunction of the product.


- Store and transport at temperatures between 0°C/32°F - 40°C/104°F
- Use the product at room temperature
- Do not freeze
- Protect from moisture

3. STEP-BY-STEP INSTRUCTIONS


3.1 Initial set-up

NOTE: Upon receiving the SmartLite® Max L.E.D. Curing Light, check the packaging and parts for any possible damage that may have occurred in transit. If damage is apparent, please contact your authorized DENTSPLY distributor from which the light was purchased.

Remove the contents from the shipping package. Be careful not to drop or impact the curing light. Choose a flat, level, sturdy surface to place the charging cradle on. The SmartLite® Max Light is shipped in a “hibernation” mode, and the display will be blank. To wake the unit, press the On/Off button. (Fig 1.)


Before using for the first time, the SmartLite® Max Light must be fully charged. For the first charge, this process will take about 3 hours. Use only the power supply furnished with the SmartLite® Max Light system. Attempting to use another power supply source may cause injury or damage and will void the SmartLite® Max Light warranty.

Attach the probe to the hand piece as shown in Figure 2.

Plug the connector of the power supply into the back of the charging cradle, or plug it into the bottom end of the SmartLite® Max Light itself, as shown in Figure 3. Plug the other end of the power supply into the appropriate electrical outlet. As mentioned above, wait for at least 3 hours to completely charge the light before the first use.

NOTE: Verify the line plug of the power supply plug is properly rated for your facility, and that all connections comply with safe practices and all local requirements. Always use unit in conjunction with a properly grounded outlet.

3.2 Built-in radiometer

The charging cradle incorporates a built-in power meter, allowing you to see how well the SmartLite® Max Light is performing. (Fig. 4)

The power meter is not intended to measure the performance of the PULSE mode. Select one of the other modes, and set the timer for 10 seconds. Place the head of the SmartLite® Max Light over the sensor and activate the light. The power meter will illuminate, showing the relative power output. In STANDARD and RAMP modes, two green bars above the green ones, indicating the BOOST mode is working at a higher power than the other modes. Yellow, red or orange lights indicate insufficient power output for curing (see Troubleshooting section).


Fig. 4

Inadequate polymerization due to insufficient curing

- Check compatibility of curing light
- Check curing cycle
- Check curing output before each procedure

3.3 Modes

There are 4 different modes to select from based on the procedure and the requirement. (Fig. 5)

- **Standard** – The light output is at a steady 1200mW/cm².

- **Boost** – This mode provides additional power to the LEDs, helping them generate a higher power output up to 1400mW/cm². Since the LEDs function at such high power level, only the 5 second curing cycle is provided to avoid premature LED failure. You can use multiple boost cycles in succession to obtain longer Boost curing.

- **Pulse** – The light output is pulsed 10 times a second, giving a net effect of half the power output as the continuous mode.

- **Ramp** - There is a steady increase of the output from 0 to 1200mW/cm² within 4 seconds and then remains constant.


Fig. 5


Timer selection

Curing time can be selected to run at 5, 10, 15 or 20 seconds for all the modes except the Boost mode which is programmed to run only at 5 seconds.

3.4 Keypad

There are 3 buttons on the key pad for the timer selection, mode selection and the on/off. (Fig. 6) The on/off button is used to activate and deactivate the LEDs. Push once to turn on the light and push it again if you need to turn it off at any time during the curing cycle. The light will default back to the full time selected if interrupted.

Fig. 6

The probe housing the two LEDs connects to the hand piece on top. The two LEDs are confined within the reflector surface. The probe also houses a heat sink that helps dissipate the heat generated by the LEDs to avoid the build up of heat that could damage the LEDs.

3.5 Battery status and cord indicator

The battery indicator on the SmartLite® Max Light shows the current status of the battery. When the RED BATTERY LOW INDICATOR blinks, this means the battery has a low charge and urgent action should be taken to either place the unit back in the charging cradle, or attach the power cord directly to the handle. The unit may still operate, but the remaining runtime is very close to 0. If this symbol is continuously on, the battery doesn't have enough power to run the light and the unit will not function.

The SmartLite® Max Light is designed so that even with a low battery, the high output power indicated in the specifications is maintained. Once the output power can no longer be sustained at that level, the SmartLite® Max Light will shut down. This way, you are assured that the SmartLite® Max Light always delivers the top light output, regardless of the battery's condition.

To charge the battery, place the SmartLite® Max Light in the charging cradled, with the power supply attached to the cradle, or attach the power supply directly to the handle. When this is done, the AMBER BATTERY CHARGING INDICATOR will illuminate.


The SmartLite® Max Light may be used with the power supply attached to the handpiece handle, even if the battery is low. While in this configuration, the battery will charge while the unit is not being used. While charging, the GREEN BATTERY STATUS INDICATOR will illuminate. When the battery is fully charged, the green indicator will switch off. When the handle is removed from the charger, or the power disconnected, none of the battery indicators will illuminate.

3.6 Operation

The SmartLite® Max Light features a “sleep” mode when not in the charging cradle or attached to the power supply. If the unit has not been used for more than 10 minutes, the display will turn off to conserve power. If the display is not illuminated, simply pick up or move the unit to bring it out of the sleep mode.

The curing wand has an arrow embossed on either side of the LED head for alignment with the target (Fig. 7). Press the MODE and TIMER buttons to select the desired settings. Since the probe can be rotated 360°, it could be used both as a wand style with thumb operation or gun style with index finger operation.

Wand Style


Fig. 7

Gun Style

Always use the settings recommended by the manufacturer of the dental material when selecting the settings for the SmartLite® Max Light.

To position the head to the desired angle, simply rotate the head (plastic) portion of the SmartLite® Max Light. The head can rotate an unlimited number of times in any direction. This helps the operator to achieve easy access to remote areas of the oral cavity without causing unwanted discomfort or stress to the operator or the patient.

Probe lens damage by contamination

- Prior to use, put a disposable barrier sleeve onto the probe and make sure it fits properly without any wrinkles.
- Do not allow unprotected lens contact with restorative material.

Make sure everyone is wearing the appropriate eye protection. Once the settings are correct, position the head of the SmartLite® Max Light over the site to be illuminated and press the START button. Once the curing cycle is complete, the unit will automatically turn the light off. If needed, the START button can be pressed again to deactivate the light before the curing cycle is complete.

Be aware that the SmartLite® Max Light does generate some heat when in use. If multiple curing cycles in rapid succession occur, especially in Boost mode, the built-in temperature overload sensor may activate, causing the Boost symbol on the SmartLite® Max Light to flash. If this occurs, place the unit in the charging cradle or on the counter and wait at least 6 minutes before trying to use the unit again. This protects the lifespan of the LEDs in the unit.

4. CARE AND MAINTENANCE

4.1 Hygiene

The SmartLite® Max Light is a state-of-the-art electronic device. To prevent unnecessary damage to the unit and to the power supply, these care and cleaning guidelines must be followed.

- Remove and discard barrier sleeve.

Cross-contamination

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Reprocess reusable products as described below.

- Do not remove curing probe from handpiece until handpiece cleaning and disinfection is completed as outlined.

4.1.1 Cleaning

- Disconnect the power supply plug from the power outlet and charging base unit prior to cleaning.
- Clean the unit by gently scrubbing with a disposable towel soaked with hot water and soap or detergent. Do not allow solution to penetrate the casing.
- Always protect lens of the light probes from abrasive materials. Should the surface become marred or scratched, light transmission and curing performance will be diminished.

4.1.2 Disinfection

- Probe, handpiece, light shield and/or charging base exposed to spatter or spray of body fluids, or touched by contaminated hands or oral tissues, should be disinfected with a hospital-level disinfectant according to national/local regulations. Water-based phenolics and dual or synergized quaternary ammoniums are approved disinfectants.
- Do not immerse units in disinfectant. Disinfect the units by spraying disinfectant solution onto a paper towel and wiping or wiping with impregnated wipes. Do not spray disinfectant directly onto the units or allow solution to penetrate the casing.
- Use of iodophor, chlorine or glutaraldehyde based disinfectant solutions will void the SmartLite® Max Light warranty.
- DO NOT immerse unit in water or disinfectant. Wipe – do not spray – solution (with curing probe in place) onto plastic parts. Prevent liquids from entering openings on unit.
- DO NOT autoclave or otherwise sterilize these components.
- Reconnect the power supply plug to the power outlet and the charging base unit when cleaning and disinfection process is completed.

4.2 Head replacement

The head portion of the SmartLite® Max Light may be removed for replacement or cleaning. To do this, simply grasp the head firmly and pull it away from the lower handle portion. Attach the replacement head by inserting the connector into the top of the handle. Make sure you push the head all the way against the handle to ensure a secure connection. Make sure the power supply is disconnected from the unit before replacing the head.

5. TROUBLESHOOTING

PROBLEM: There is no display.

CAUSE: The unit is in hibernation mode.

SOLUTION: Press the on/off button to wake the unit.

PROBLEM: The unit does not appear to have any power.

CAUSE: The battery is drained.

SOLUTION: Attach the power supply to the SmartLite® Max Light handle, or place the SmartLite® Max Light in the charging cradle, making sure the power supply is attached to the cradle.

CAUSE: The power supply is not attached to the SmartLite® Max Light or the electrical outlet.

SOLUTION: Make sure the power supply is securely connected to electrical supply outlet and to the SmartLite® Max Light handle or cradle.

CAUSE: The SmartLite® Max Light is in sleep mode.

SOLUTION: Press any button once to bring the unit out of sleep mode. The display should illuminate when this is done.

PROBLEM: The head or handle is getting hot.

CAUSE: Unit has been used for too many consecutive cycles.

SOLUTION: Place the SmartLite® Max Light in the charger, or set the unit on a counter, and allow it to cool down for at least 6 minutes.

PROBLEM: Power output indicated on the cradle does not match the results from another power meter.

CAUSE: The cradle power meter is not a quantitative measuring device.

SOLUTION: The power meter incorporated into the charging cradle is meant as a simple diagnostic tool only, and is calibrated to the output performance of the SmartLite® Max Light. Other curing lights may not give a similar result on the SmartLite® Max Light power meter as the SmartLite® Max Light does. Due to the oval shape and dispersion angle of the SmartLite® Max Light, the power measured on another power meter may not be accurate.

PROBLEM: The power meter doesn't show a green bar when the output is measured.

CAUSE: SmartLite® Max Light head is not placed in the sensor correctly.

SOLUTION: Make sure the SmartLite® Max Light head is resting securely in the sensor, with the head resting against the sensor face.

CAUSE: SmartLite® Max Light is in Pulse mode.

SOLUTION: Do not measure the power output in Pulse mode. Select another mode to measure the power output.

PROBLEM: I can't select the time in Boost mode.

CAUSE: Timer function is not available in Boost mode.

SOLUTION: If more time is desired, either use the Boost mode repeatedly at the 5 second curing cycle, or select another function and choose the desired curing time.

PROBLEM: I press Start, but no light comes out.

CAUSE: Head is not attached properly.

SOLUTION: Make sure the head is pushed completely into the handle, such that the head connector is not visible.

PROBLEM: I've been using the unit for several cycles and it suddenly stopped working.

CAUSE: Temperature overload – sensor has activated.

SOLUTION: The SmartLite® Max Light is equipped with a temperature overload sensor to protect the LEDs from excessive heat. If this occurs, place the unit in the cradle or on the counter and wait at least 6 minutes for the unit to cool.

CAUSE: Battery has been drained.

SOLUTION: The battery has been drained of its stored energy. Place the unit in the charging cradle, or attach the power supply to the handle.

PROBLEM: The red battery low indicator light is on and I plugged in the power cord, but the light still won't function.

CAUSE: The light has to reset when the battery low indicator light is on continuously and the power cord is plugged in for corded use.

SOLUTION: plug in the power cord and let the light reset for several minutes before use.

If the above steps do not solve the problem, contact the dealer or point of purchase. DO NOT USE THE SmartLite® Max Light AND DO NOT ATTEMPT TO SERVICE THE UNIT. Serious injury may result.

SmartLite® Max L.E.D. Curing Light complies with the following Directives/Standards:

- Medical electrical equipment Part 1: General requirements for safety (IEC 60601-1)
- Standard for Safety for Medical Electrical Equipment, Part 1: General Requirements for Safety (UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90)
- Medical Electrical Equipment - Part 1-4: General Requirements for Safety - Collateral Standard: Programmable Electrical Medical Systems Consolidated (IEC 60601-1-4)

6. WARRANTY

This product is designed for use in a dental office and this warranty is not applicable to other uses. DENTSPLY warrants this product against defective materials and workmanship for a period of twenty four (24) months from the date of purchase. In the event of such a defect, DENTSPLY will repair or replace the product or necessary parts therein, at its discretion, and such repair or replacement shall be the sole remedy of this warranty. This warranty extends only to the original purchase and is subject to these conditions:

- The unit must not be subjected to abuse, improper installation or application, nor repair service by other than trained DENTSPLY service personnel.

THERE ARE NO WARRANTIES, EXPRESS OR LIMITED, WHICH EXTEND BEYOND THIS DESCRIPTION. DENTSPLY neither assumes, nor authorizes any person to assume for it, any other liability in connection with the sale and use of this product.

DAMAGES ARE LIMITED STRICTLY TO REPLACEMENT OF THE PRODUCT. DENTSPLY EXPRESSLY DISCLAIMS LIABILITY FOR INCIDENTAL AND CONSEQUENTIAL DAMAGES RESULTING FROM ITS USE.

Claims covered by this warranty will be honored when this warranty is presented within 30 days from discovery of defect. During the warranty period, DENTSPLY International reserves the right to repair or replace this product.

7. SERVICE

For product service please contact DENTSPLY Caulk customer service at 1-800-532-2855. Customers outside the USA should contact their nearest DENTSPLY location, as listed on the back cover of this manual.

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DENTSPLY Limited
Aviator Park
Station Road
Addleston, Weybridge
Surrey KT15 2PG
Tel. 44 01932 853422

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DENTSPLY Italia S.r.l.
Curtatone, 3
I-00185 Rome
Italy
Tel. 39 06 72 64 031

DENTSPLY Australia Pty. Ltd.
11-21 Gilby Road
Mount Waverley Victoria 3149
Australia
Tel. 61-03-9538-8280

DENTSPLY Sankin
Azabu Kaisei Building
1-8-10, Azabudai, Minato-ku
Tokyo, 106-0041
Japan
Tel. 81-03-514-1001

Hong Kong Division of DENTSPLY International
25/F, Pacific Link Tower
Southmark Tower
11 Yip Hing Street
Wong Chuk Hang
Hong Kong
Tel. (852)-2870-0336

DENTSPLY Industria e Comercio Ltda.
Praia de Botafogo
228/Bloco A/Grupo 1501-B
Botafogo
22359-900 Rio de Janeiro, RJ
Brazil
Tel. (55-21)-3237-9800

DENTSPLY Mexico, S.A. de C.V.
Calzada Vallejo No. 846
Col. Industrial Vallejo
02300 Mexico, D.F.
Tel. 52-55-5587-6488

DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3 Canada
Tel. 800-263-1437

Manufactured for:
DENTSPLY Caulk
38 West Clarke Avenue
Milford, DE 19963 USA
1-302-422-4511
www.dentsply.com

Distributed by:
DENTSPLY Canada
161 Vinyl Court
Woodbridge, Ontario
L4L 4A3 Canada